



AVENUE SYNC

Easy solution to support Aurora IQ handheld fundus camera integration to PACS and worklist server.

OPTOMED CORPORATION

Avenue Sync features




Avenue Sync integrates Aurora IQ camera seamlessly with hospital systems and enables more efficient work with existing tools.



Avenue Sync transfers:

- DICOM worklist from hospital information system to Aurora IQ camera
- DICOM images from Aurora IQ camera to hospital image archive.



Worklist

Image



DICOM Worklist

DICOM Image

Extra features



Tools for adding, editing and deleting study information



Supports different workflows



Possibility to view images



Simple UI for effective use